

600 mm Helmholtz Coils

Ferronato® - BH600-3-A

- --- Desktop set of three pairs of Helmholtz coils, for laboratory and general purposes.
- --- Suited for many magnetic measurements and experiments, in DC and AC.
 - Equal generating field rate for the three pairs: 216 μT/A
 - Accurately made, with error smaller than ±1% in the generated field.
 - Thanks to its simple support and its joins by screws the coil arrengement can be modified with relative facility.
 - · Coils on aluminium alloy forms.
 - Each aluminium form provides a usable extra turn, with connection in the terminal block. An application example is the generation of a small magnetic field (DC or AC) to modulate the main one. Also it can be wired to generate small gradients.
 - The aluminium forms also act like electrostatic screens.
 - The coils can undergo heatings of until at least 100 °C without damage.
 - Robust construction but with a reasonable weight.
 - Totally constructed with non-ferromagnetic materials.
 - Excellent quality/price ratio.
 - There are versions of one and two axis available, with similar characteristics:
 - o BH600-2-A Set of two pairs (two axes), X and Y.
 - o BH600-1-A Set of one pair (one axis), X.







SPECIFICATIONS OF THE SET (BH600-3-A)

Field/Current ratio:	216 μT/A (2,16 Gauss/A). For each pair, X, Y or Z. Maximum error: ±1%.		
Maximum field:	864 μT (8.64 Gauss) in permanent mode / 2.2 mT (22 Gauss) during 2 minutes. Each pair.		
Maximum current:	4.0 A in permanent mode / 10 A during 2 minutes (start temperature: 20 °C). Each pair.		
Isolation voltage:	250 V DC minimum, between windings and their forms and between pairs. Tested to 500 V DC.		
Magnetic field homogeneity:	Differences smaller than ±1% with respect to the centre, in a spherical volume of 150 mm of diameter centred in the coils. Differences smaller than ±5% in a spherical volume of 220 mm of diameter. Volumes to 1% and 5% greater in some directions.		
Connection:	Connection block with Ø4 mm screws (M4).		
Maximum working temperature:	80 °C for the set / 100 °C for the coils.		
Coil cross section:	Winding: 8.5 x 10 mm, maximum. Total (form): 10 x 13 mm		
Materials:	Copper winding. Coil forms of aluminium alloy, with connection blocks of resin/glass fibre (FR4) with PVC covers. Stand support of PVC (rigid PVC in pillars, foamed PVC in plates superior and inferior), with brackets of Acetal. Screws of nickel-plated brass and Nylon.		
Maximum dimensions:	Height 678 mm x Wide 609 mm x Depth 576 mm.		
Weight:	10.8 kg for the BH600-3-A set (smaller for BH600-1-A/2-A).		
Accessories:	Delivered with Instructions of Use in Spanish and English.		
Warranty:	Two years.		

SPECIFICATIONS FOR EACH COIL PAIR

	X pair (large)	Y pair (medium) (2)	Z pair (small) (3)
Effective diameter:	600 ±1 mm	567 ±1 mm	533 ±1 mm
Turns number:	72	68	64
DC Resistance, at 20 °C: (1)	5.9 Ω ±3%	5.3 Ω ±3%	4.7 Ω ±3%
Self-inductance:	18.2 mH ±5%	15.2 mH ±5%	12.6 mH ±5%
Self-resonance Frequency: (With one form end connected to one coil end)	~ 16 kHz	~ 17 kHz	~ 19 kHz
Secondary field generated by the forms when used as coils (Xs, Ys, Zs): (4)	2.99 μT/A ±1%	3.17 μT/A ±1%	3.37 μT/A ±1%

^{(1) -} Resistance measured at the general connection block.

For any enquiry, please, do not hesitate in contact us:

Our Telephone: (+34) 925 536154Our Fax: (+34) 925 537644

• Our E-mail: serviciencia@serviciencia.es

• Internet: www.serviciencia.es

^{(2) -} Pair not included in the BH600-1-A set version.

 $^{(\}ensuremath{^3})$ - Pair not incluided in the BH600-1-A and BH600-2-A set versions.

^{(4) -} We call this constuctive idea "In-Circuit Coils Forms".

⁻ These specifications are subject to change without prior notice -

⁻ The set is supplied mounted and ready to use, in a wooden box.