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
## 2.- FIGURE 1 EXPLANATION

- 1) Green pilot. Blinks on every QRS peak. Indicates power-on and rate.
- 2) Push-button to switch 60 to 130 bpm and vice versa.
- 3) Red pilot. Starts to blink when the battery should be replaced.
- 4) Push-button for simulator switch-on.
- 5) Push-button for simulator switch-off.
- 6) Battery lodging.
- 7) Output connectors. Colours according to IEC.
- 8) Detail shows how to push the cover tab, to open it and access the battery.

## 3.- FC12D TECHNICAL SPECIFICATIONS.

<b>Waveform:</b>	"Normal" PQRST (see Fig. 1 and 2 on page 2).
<b>Output Amplitude:</b>	Fixed, 1.0 mV peak-peak ( $\pm 5\%$ ), for Lead I, at 60 bpm (see Paragraph 4 for the others).
<b>Signal Rate:</b>	Two selectable values: 60 or 130 bpm ( $\pm 0.2\%$ ). (bpm = beats per minute).
<b>Output Impedance:</b>	4500 $\Omega$ , maximum, between any terminal and N (RL).
<b>Output Connectors:</b>	Ten of 4-mm "Banana" type and ten of "Snap-On" type, in parallel.
<b>Indicating Lights:</b>	Two. The Green one blinks on each QRS. The Red one blinks under "Low Battery" condition.
<b>Power Supply:</b>	9 V battery, PP3 or 6LR61 (alkaline) type.
<b>Power Consumption:</b>	310 $\mu\text{A}$ (integrated value), at 60 bpm.
<b>Operating Temp.:</b>	From $-10^{\circ}\text{C}$ to $+50^{\circ}\text{C}$ .
<b>Size:</b>	Maximum: width 74 x length 138 x depth 31 mm.
<b>Weight:</b>	140 grams without battery. About 185 g with battery.

## 4.- OUTPUT AMPLITUDES (in mV peak to peak, $\pm 5\%$ )

I = 1.00	aVR = 1.12
II = 1.25	aVL = 0.37
III = 0.25	aVF = 0.75

The Leads V1 to V6 could have a bigger amplitude variation from one simulator to another, estimated between +10% and -20% in respective to shown in Fig. 2.

## 5.- THE BATTERY.

5.1- Access: to install it, or to replace it, it is necessary to remove its lodging cover, with a tool to push the tab, as seen in Fig. 1. So, the battery and its connector are attainable.

5.2- Battery type: it is a 9 Volts one. It is possible to use the lower priced PP3 type, but it is advisable to use alkaline ones (6LR61 type), for a much longer life and a lower risk of chemical leaks.

5.3- Battery life: it depends on the employed type. The alkaline ones can yield several hundreds of working hours (300 or more).

5.4- "Low Battery" indication: the Red pilot (3 in the drawing) starts to blink when battery voltage goes down to about 6.3 V. At that moment it is advisable to replace the battery for a new one. The FC12D output signal is sustained without much amplitude variations until the battery voltage reaches about 6.0 V, under which the simulator will work for some more time but signal amplitude will decrease.

## 6.- THE FC12D USE.

6.1- Switch-On: a delay was introduced to reduce the risk of involuntary actions. It is necessary to keep the "ON" push-button (5 in drawing) pressed for about 3 seconds, or until the "QRS" pilot (1) starts to blink. This pilot lights in coincidence with the QRS peak. The simulator will always be at 60 bpm after switch-on.

6.2- Switch-Off: "OFF" (4) must be pushed.

6.3- 60 to 130 bpm switching: to push on "60/130 bpm" (2) for about 3 seconds. The same operation must be done to switch from 130 to 60 bpm. The "QRS" pilot keeps dark while the push-button is pressed.

6.4- Connection: to get the specified signal amplitudes the cables from the equipment under test must be properly attached to the simulator. Otherwise the signal could be inverted, or have double amplitude, or be noisy, etc.

## 7.- THE FC12D OPERATION.

7.1- Waveform generation: it is performed by one specially programmed microcontroller ( $\mu\text{C}$ ), in which memory are the values to shape the main wave. The oscillator on the  $\mu\text{C}$  is controlled by one quartz crystal, due to which the signal rhythm is very stable. The  $\mu\text{C}$  also receives the commands from the push-buttons and performs the corresponding action. The digital output from the  $\mu\text{C}$  is transformed into an analogue one by one D/A converter circuit and it is then slightly filtered. Another circuit generates the big "S" peak negative, to be added to the main complex in order to get the precordial signals (V1...V6).

7.2- Power supply: the battery voltage is regulated by one specialised microcircuit to ensure one constant signal amplitude not depending on the battery voltage, which also performs as a threshold detector to indicate the "low battery" situation.

7.3- Adjustments: there is no prevision for adjustments after fabrication. However, if it is necessary, there is one signal amplitude adjustment available, to be done by the pre-set potentiometer visible on the printed circuit when the back cover is retired.

7.4- Waveforms at 60 or 130 bpm: these are not exactly the same. Some internal periods in the PQRST complex are shortened at 130 bpm, as well as its total time width, in the way as it happens in the human body. That can be seen in Fig. 3.

## 8.- GENERAL CONSIDERATIONS ABOUT THE FC12D USE

8.1- The simulator has ten duplicated output connections, each one with two types of connectors: 4-mm female Banana and Snap-On plug. The Bananas are normally found on electrocardiograph connection cables, meanwhile Snap-On are on monitors, Holter recorders and equipment using ECG as trigger signal for other functions.

8.2- The connector named N (or RL) might be, at first, be always connected to the same name terminal on the equipment under test, or to the Common (or "Neutral") in case of one Holter recorder.

8.3- Following are explained the standard terminal names and the patient area in where electrocardiograph leads connections are done.

**N** = RL = Right Leg (right ankle) = "Neutral" = Reference.

**R** = RA = Right Arm (right wrist).

**L** = LA = Left Arm (left wrist).

**F** = LL = Left Leg (left ankle).

**C1 to C6** = V1 to V6 = electrodes in precordial area of chest, giving the Unipolar Chest Leads (V1 to V6).

(Due to practical reasons the electrodes of ECG monitors and other equipment are all connected to the human body trunk, instead of the extremities).

8.4- For terminal colours the international code for electrocardiographs, from the International Electrotechnical Commission (IEC) is followed (sometimes termed as "European"), but it is possible to find other colours on equipments. On little known equipment it is advisable to do the connections according to its terminal names, instead of its colours. In case of doubt some experiments should be convenient until getting the expected signal. That procedure can be useful to identify non marked terminals. The most important terminal to be identified on the equipment under test is N, or RL, since it is the reference for all the signals.

## **9.- USUAL APPLICATIONS**

### **9.1.- With cardiac monitors or electrocardiographs.**

9.1.1- This is the easiest case. Simply, the terminals are connected to the respective FC12D connectors.

9.1.2- The signal shown on the monitor (or electrocardiograph) might be of 1,0 mVp-p for Lead I, or 1,25 mVp-p for Lead II, etc., with a rate of 60 or 130 bpm according the selected on the simulator.

### **9.2.- With one channel Holter recorders (three terminals).**

9.2.1- The Neutral terminal, usually of Green colour, is attached to N (RL) on the simulator. The non-inverting (+) terminal, usually Red, is attached to L (LA) on the simulator and the inverting terminal (-), usually White, to R (RA).

9.2.2.- With the above connections the signal will be 1,0 mVp-p. The signal is inverted when the (+) and (-) terminals connection is inverted.

### **9.3.- With two channels Holter recorders (five terminals).**

9.3.1- The procedure for Channel 1 is same as indicated in 9.2.1.

9.3.2- For Channel 2 the terminals could be hooked-up, by means of one adapter, together with those from Channel 1. So, the (+), usually Black to connector L and the (-), usually Brown, to connector R. The signal will be the same for the two channels.

9.3.3- One option to connect the Channel 2 without one adapter is to use the precordial outputs. So, if the (+) (Black) is attached to C5 and the (-) (Brown) to C1, a signal is provided with the same shape as in Channel 1 and with a similar amplitude (a little wider).

**9.4.- With three channels Holter recorders (seven terminals).**

9.4.1- Channels 1 and 2 can be connected as said in 9.3.

9.4.2- For Channel 3 the terminals could be hooked-up, by means of one adapter, together with those from Channels 1 and 2. So, the (+), usually Orange to connector L and the (-), usually Blue, to connector R. The signal will be the same for the three channels.

9.4.3- Without one adapter and with Channel 2 connected according 9.3.3, it is possible to hook-up Channel 3 to other precordial outputs. So, if (+) (Orange) is attached to C6 and (-) (Blue) to C2, a signal is provided which is very similar to Channel 2 one.

**9.5.- With other various equipment.**

9.5.1- There are much equipment using ECG as a trigger signal. Commonly a 3 terminal connection is used and sometimes a 4 or 5 terminal one. The previously mentioned for monitors (9.1) is applicable here. Some tests with different ways in hooking-up would allow the configuration for a suitable signal to be found in case of doubt.

**DECLARATION OF CONFORMITY**

We, SERVICIENCIA, S.L. declare under our sole responsibility that this apparatus conforms to the Council Directive **89/336/EEC**, amended by **91/263/EEC** and **92/31/EEC**, about EMC. In particular it conforms to the following harmonised EMC standard:

EN 61326-1:1997+A1(1998); Emissions and Immunity for Equipment for Measurement, Control and Laboratory Use.

Alberto R. Marino  
General Manager  
Serviciencia, S.L.  
07 January 2002

Manufacturer: SERVICIENCIA, S.L.

Isabel II, 22 - 45210 YUNCOS (Toledo) - SPAIN

Phone: (+34) 925 536154 – Fax: (+34) 925 537644 - Internet: [www.serviciencia.es](http://www.serviciencia.es)