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

1.- FC3D TECHNICAL SPECIFICATIONS.

Waveform:	"Normal" PQRST.
Output Amplitude:	Fixed, 1.0 mV peak-peak ($\pm 5\%$), for Lead I, II or III, at 60 bpm.
Signal Rate:	Two selectable values: 60 or 130 bpm ($\pm 0,2\%$). (bpm = beats per minute).
Output Impedance:	4000 Ω , maximum.
Output Connectors:	Three of 4-mm Banana type, female, 20 mm spaced (colours according to Lead I) and three of "Snap-On" type studs, in parallel.
Indicating Lights:	Two. One (Green) blinks on each QRS. Another (Red) blinks under "Low Battery" condition.
Power Supply:	9 V battery, PP3 or 6LR61 (alkaline) type.
Power Consumption:	310 μA (integrated value), at 60 bpm.
Operating Temp.:	From -10°C to 50°C .
Size:	Maximum: width 64 x length 91 x depth 26 millimetres.
Weight:	65 grams without battery. Approximately 110 g with battery.

2.- WAVEFORMS.

Images obtained with a 0 to 370 Hz (-3db) bandwidth ECG amplifier.

3.- THE BATTERY.

3.1- Access: to install or to replace it, it is necessary to remove the screw on the backside and to lift the cover with the help of the notch on the screw side. So, the battery and its connector are attainable.

3.2- Battery type: it is a 9 Volt one. It is possible to use the lower priced PP3 type, but it is advisable to use alkaline ones (6LR61 type), for a much longer life and a lower risk of chemical leaks.

3.3- Battery life: it depends on the employed type. The alkaline type can yield several hundreds of working hours (300 or more). It is advisable to remove the battery when the device will be unused for a long time.

3.4- "Low Battery" indication: the Red pilot (3 in the drawing) starts to blink when its voltage goes down to about 6.3 V. At that moment it is advisable to replace the battery for a new one. The FC3D output signal is sustained without much amplitude variations until the battery voltage reaches about 6 V, below which the simulator will work for some more time but signal amplitude will decrease.

4.- FC3D USE.

4.1- Switching-On: a time delay was introduced to reduce the risk of involuntary actions. It is necessary to keep the "ON" push-button (5 in drawing) pressed for about 3 seconds, or until the "QRS" pilot (1) starts to blink. This pilot lights in coincidence with the QRS peak. The simulator will always be at 60 bpm after switch-on.

4.2- Switching-Off: "OFF" (4) must be pushed.

4.3- 60 to 130 bpm switching: push on "60/130 bpm" (2) for about 3 seconds. The same operation must be done to switch from 130 to 60 bpm. The "QRS" pilot stays dark while the push-button is pressed.

4.4- Connection: to get the specified signal amplitude (1.0 mVp-p) the cables from the equipment under test must be properly attached to the simulator, according to what is shown on the FC3D front face for each Lead (I, II or III). Otherwise the signal could be inverted, or to have double amplitude, or to be noisy, etc.

5.- FC3D OPERATION.

5.1- Waveform generation: it is performed by a specially programmed microcontroller (μC), in which memory are the values to shape the wave. The oscillator on the μC is controlled by one quartz crystal, due to which the signal rhythm is very stable. The μC also receives the commands from the push-buttons and performs the corresponding actions. The digital output from the μC is transformed into an analogue signal by a D/A converter circuit and, then, it is slightly filtered.

5.2- Power supply: the battery voltage is regulated by one specialised microcircuit to ensure one constant signal amplitude not depending on the battery voltage, which also performs as a threshold detector to indicate the "low battery" situation.

5.3- Adjustments: there is no prevision for adjustments after fabrication. However, if it is necessary, there is one signal amplitude adjustment available, to be done by the pre-set potentiometer visible on the printed circuit when the back cover is removed.

5.4- Waveforms at 60 or 130 bpm: these are not exactly the same. Some internal periods in the PQRST complex are shortened at 130 bpm, as well as its total time width, in the way as it happens in the human body.

6.- GENERAL CONSIDERATIONS ABOUT THE FC3D USE.

6.1- The simulator has three duplicated output connections, each one with two types of connectors: 4-mm female Banana and Snap-On stud. The bananas are normally found on electrocardiograph connection cables, meanwhile Snap-On are on monitors, Holter recorders and equipment using ECG as a trigger signal for other functions.

6.2- The connector named N (or RL) might be, at first, always connected to the same terminal on the equipment under test, or to the Common (or "Neutral") in case of one Holter recorder. The other two connectors, numbered 2 and 3, accept different terminals, according to the wanted lead. On the FC3D front face are indicated the terminals to be connected in 2 and 3 to get the signal for Leads I, II and III (which are the possible ones with three terminals). By means of one adapter it is possible to connect the 10 terminals from one electrocardiograph, as explained in below (7.4).

6.3- The standard terminals denomination is explained in the following, as well as the corresponding patient areas where the electrodes from one electrocardiograph are connected.

N = RL = Right Leg (right ankle) = "Neutral" = Reference.

R = RA = Right Arm (right wrist).

L = LA = Left Arm (left wrist).

F = LL = Left Leg = (left ankle).

C1 to C6 = V1 to V6 = electrodes in precordial area of chest, giving the Unipolar Chest Leads (V1 to V6).

(Due to practical reasons the electrodes of ECG monitors and other equipment are all connected to body trunk, instead of the extremities).

6.4- For terminal colours the international code for electrocardiograph, from the International Electrotechnical Commission (IEC), sometimes termed as "European", is followed, but it is possible to find other colours on equipment. On little known equipment it is advisable to do the connections according to its terminal names, instead of its colours. In case of doubt some experiments should be convenient until getting the expected signal. The most important terminal to be identified on the equipment under test is N, or RL, since it is the reference for all the signals.

7.- USUAL APPLICATIONS

7.1- With three terminal cardiac monitors.

7.1.1- This is the easiest case. Simply the terminals are connected to the respective FC3D connectors.

7.1.2- The signal showed on the monitor might be of 1.0 mVp-p and with a rate of 60 or 130 bpm, according to the selected on the simulator.

7.1.3- Problems: in case of inverted terminals on connectors 2 and 3, the signal will also be inverted. If the monitor terminal N (RL) were in connector 2 the signal will have double to the specified amplitude (and perhaps more noise).

7.2.- With four terminal monitors

7.2.1- The most important is to locate the monitor N (RL) terminal and connect it to the simulator corresponding connector.

7.2.2- To have signal on the three Leads I, II and III it is possible to connect the terminals according to Lead I and, by means of some adapter, the spare F (LL) terminal to connector 1, together with N terminal. The signal on Lead II will be inverted and Lead III will be inverted and with double amplitude.

7.2.3- To have signal on the three Leads aVR, aVL and aVF it is possible to make the connection for Lead I (like before, on 7.2.2) and the remaining terminal to connector 3, together with L (LA) terminal.

7.3.- With a five terminal monitor.

7.3.1- To have signal on Lead V it is possible to connect the terminals as in 7.2.2. and 7.2.3., connecting the remaining terminal C in the connector 3, together with L (LA) terminal.

7.4.- With a electrocardiograph.

7.4.1- In general it is true the previously mentioned in 7.2 and 7.3, for 5 terminal monitors.

7.4.2- If the electrocardiograph has 10 terminals it is possible to get signals (all the same) for Leads V1 to V6 by connecting C1 to C6 terminals, all together or one to one alternately, on the connector 3, together with L (LA) terminal. See 7.3.1.

7.5.- With one channel Holter recorder (three terminals).

7.5.1- The Neutral terminal, usually Green, is connected to N (RL) on the simulator. The non-inverting terminal (+), usually Red, is connected to connector 3 on the simulator and the inverting terminal (-), usually White, to connector 2.

7.5.2- With the connections as above there will be a 1,0 mVp-p signal. The signal gets inverted when terminals (+) and (-) connections are inverted.

7.6.- With two channel Holter recorders (five terminals).

7.6.1- For channel 1 the procedure is the same as indicated in 7.5.1.

7.6.2- For Channel 2 the terminals can be connected, by means of one adapter, together with those from Channel 1. So the (+) usually Black, to connector 3 and the (-), usually Brown, to connector 2. The signal will be the same on both channels.

7.7.- With three channel Holter recorders (seven terminals).

7.7.1- Channels 1 and 2 can be connected according to the previously mentioned in 7.6.

7.7.2- For Channel 3 the terminal can be connected, by means of one adapter, together with those from Channels 1 and 2. So the (+), usually Orange, to connector 3 and the (-), usually Blue, to connector 2. The signal will be the same on the three channels.

7.8.- With other various equipment.

7.8.1- There are much equipment using the ECG as a trigger signal. Commonly a 3 terminal connection is used and sometimes a 4 terminal one. In case of 3 terminals these could be connected as indicated on 7.1.1, for simple monitors. In the case of 4 terminals these could be connected as said in 7.2.2. Some tests with different ways in hooking-up would allow the configuration for a suitable signal to be found in case of doubt.

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DECLARATION OF CONFORMITY

We, *SERVICIENCIA, S.L.* declare under our sole responsibility that this apparatus conforms to the Council Directive **89/336/EEC**, amended by **91/263/EEC** and **92/31/EEC**, about EMC. In particular it conforms to the following harmonised EMC standard:

EN 61326-1:1997+A1(1998); Emissions and Immunity for Equipment for Measurement, Control and Laboratory Use.

Alberto R. Marino
 Administrador / *General Manager*
 Serviciencia, S.L.
 07 Enero 2002 / *07 January 2002*

Fabricante / *Manufacturer:* SERVICIENCIA, S.L.

Dirección / *Address:* Isabel II, 22
 45210 YUNCOS
 (Toledo)
 ESPAÑA / *SPAIN*

Tel/*Phone:* (+34) 925 536154 – Fax: (+34) 925 537644 - Internet: www.serviciencia.es